






High Quality Dental Instruments



☎ 646 688 5552

🌐 cureventory.com

✉ info@cureventory.com

📍 252 GlenHill Ave, Yonkes, NY 10701

A Little About Us:

Cureventory.com Present This New Edition Dental Instruments Catalogue With Hope That It Would Be An Excellent Guide To All Over Importers, Distributors And Dental Professionals Worldwide.

Cureventory.com is ISO-9001-2008, ISO-13485-2012, CE, & cGMP (FDA-USA) Certified Company And Head Office Is Located In Industrial Area Of Sialkot, Which Is Well Known Worldwide For Manufacturing Quality Surgical & Dental Instruments. We At **Comforts** Successfully Running Family Owned Business For Last 50years With The Support Of Great Team, Whose Skilled Techniques And Traditional Approach Towards Toward Modern Surgery, Make All Possibilities To Acclaim Us Quality Instruments Manufacturer.

We Are Also Very Grateful For The Loyal Support Given By Our Customers/Distributors Who Have Been Selling Our Instruments In Their Market With Pride And Honor.

We Are Pledge To Serve More In The Years Ahead To Our Quality Minded Customers With Much More Higher Standard Of Instruments And Services.

We Offer A Comprehensive Portfolio Of Reusable And Disposable Instruments For Dental Surgeries, Dental Laboratories, And For Dental Students.

We Also Welcome Your Enquires For Instruments Not Listed In This Catalogue Both In Surgical & Dental Instruments.

At Cureventory.com We Feel That Good Communication Is Key To Any Successfully Business Relationship. Please Feel Free To Contact Us By The Methods That Are The Most Convenient For Your Business And We Will Do Our Best To Always Provide Great Quality Products And Services.

Yours Sincerely,











Cureventory.com

Cureventory.com warranties all its products against poor craftsmanship and materials. Any product found defective will be repaired or replaced free of charge. This warranty does not cover any products which become defective through misuse, negligence and mishandling.


Cureventory.com Instruments



Index



01	Extraction Forceps(English American Children).....		Page 01-09
02	Bone Rongeurs		Page 10
03	Root Elevators		Page 11-16
04	Dental Anesthesia Syringes		Page 17
05	Handles for Mouth Mirrors		Page 18
06	Dental Tweezers		Page 19
07	Dental Scissors		Page 20
08	Needle Holders		Page 21
09	Explorers & Probes		Page 23
10	Excavators / Pluggers / Condensers / Carvers		Page 24
11	Burnishers / Hollenback / Carving Knives		Page 25
12	Restorative		Page 26-27
13	Dental Scalars & Curettes		Page 28-30
14	Periodontal Knives		Page 31
15	Dental General Instruments		Page 32
16	Rubber Dam Instruments		Page 33
17	Dental Crown Instruments		Page 34
18	Dental Implants Instruments		Page 35-37
19	Orthodontics Instruments		Page 38-40
20	Impression Trays		Page 41
21	Dental Hollow Wares / Cassettes		Page 42-43





Requirement of the Dental Forceps INTRA-ALVEOLAR EXTRACTION

Dental extraction forceps has three components which is blades and handles united together by a hinge joints.
1-The serrated handle must possess a suitable size to rest comfortably in the operator's hand and should long enough to afford use to strong and steady extraction movements. Also it must be suitable in shape and designs to suit the area of the tooth.




Extracting Forceps (Set of 6 Pcs)

Universal Atraumatic Extraction Forceps Set Premium
(Adult / Pediatric / Children)

Set of 6

Cross serrations for a non-slip grasp. Easily enter the sub gingival.


Unique Design.



CSI-001



CSI-002 Set of 6 Pcs



Premium Natural Atraumatic Universal Extraction Forcep Lower Anterior
Premium Natural Atraumatic Universal Extraction Forcep Upper Anterior
Premium Natural Atraumatic Universal Extraction Forcep Universal Upper Premolars
Premium Natural Atraumatic Universal Extraction Forcep Universal Lower Premolars
Premium Natural Atraumatic Universal Extraction Forcep Lower Molars
Premium Natural Atraumatic Universal Extraction Forcep Upper Molars

Extracting Forceps (English Pattern)

Oral Surgery

Extracting Forceps (Set of 6 Pcs)

Cross serrations for a non-slip grasp. Easily enter the sub gingival.

Universal Atraumatic Extraction Forceps Set

Unique Design.



Extracting Forceps (Set of 6 Pcs)



Cross serrations for a non-slip grasp. Easily enter the sub gingival.

Black Anti-Reflection Coated

Universal Atraumatic Extraction Forceps Set

Unique Design.





CSI-101
Fig.1
Upper Centrals
and Canines




CSI-102
Fig.2
Upper Laterals
and Canines



CSI-103
Fig.3
Upper Incisors
and Canines



CSI-104
Fig.4
Lower Incisors
and Canines



CSI-105
Fig.7
Upper Premolars




CSI-106
Fig.8
Lower Premolars




CSI-107
Fig.13
Lower Premolars



CSI-108
Fig.17
Upper Molars
Right



CSI-109
Fig.18
Upper Molars
Left



CSI-110
Fig.19
Upper Third
Molars



CSI-111
Fig.20
Lower Third
Molars




CSI-112
Fig.21
Lower Molars



CSI-113
Fig.22
Lower Molars



CSI-114
Fig.30
Upper Roots



03




CSI-115
33-Lower Roots




CSI-116
33A-Lower Roots




CSI-117
73A-Lower Molars



CSI-118
79-Lower Third Molars



CSI-119
45-Lower roots



CSI-120
44-Upper roots



CSI-121
51A-Upper Roots



CSI-122
67A-Upper Third Molars



CSI-123
74-Lower Roots



CSI-124
74N-Lower Roots




CSI-125
Routurier (# 22 1/2 L)
Lower Molars and
wisdoms, left

**WISDOM
EXTRACTING
FORCEPS**



CSI-126
Routurier (# 22 1/2 R)
Lower Molars and
wisdoms, right




CSI-127
#73-Lower Molars



CSI-128
#86A-Lower Molars

DEEP GRIPPING FORCEPS

- Parallel gripping beak
- Aggressive edge serrations
- Long axis serrations within the beaks
- Thick beaks for superior penetration
- Contact the root remaining coronal structure
- Used on partially intact, broken down or sectioned teeth



No Fracture | No Tissue Damage | Extraction with Easy | Slim Light Weight Pattern



**CSI-133
35**
Upper Premolars
Deep-gripping



Thin Tip

**CSI-129
34**
Upper Incisors and
Premolars
Deep-gripping




Wide Tip

**CSI-130
34A**
Upper Centrals / Cusps
Deep-gripping




Wide Tip

**CSI-131
35A**
Upper Incisors and
Premolars
Deep-gripping



Thin Tip

**CSI-132
36**
Lower Premolars
Deep-gripping



Wide Tip

**CSI-134
36A**
Lower Premolars
Deep-gripping



**CSI-135
336**
Lower Premolars
Deep-gripping



**CSI-136
F-7C**
Lower Molar




**CSI-137
9**
Apex Forceps
Lower Universal
Deep Root Pick




**CSI-138
79N**
Lower Third Molars
Deep-gripping




**CSI-139
79N**
Lower Third Molars



**CSI-140
67**
Upper Third Molars



**CSI-141
67A**
Upper Third Molars




**CSI-142
67N**
Upper Third Molars
Deep-gripping

Extracting Forceps (English Pattern)

Oral Surgery


Extracting Forceps (Set of 10 Pcs)



05



CSI-144 Fig.1, Upper Laterals and Canines




CSI-144 Fig. 7 Upper Premolars



CSI-144 Fig. 13 lower Premolars




CSI-144 Fig. 17 Upper Molars Right



CSI-144 Fig. 18 Upper Molars Left




CSI-144 Fig. 22 Lower Roots



CSI-145 Fig. 51A Lower Roots



CSI-146 Fig. 67 Lower Third Molars



CSI-147 Fig. 74N Upper Molar




CSI-148 Fig. 79 Upper Third Molars



CSI-149 Set of 10 Pieces Fig. 1, Fig. 7, Fig. 13, Fig. 17 Fig. 18, Fig. 22, Fig. 51A, Fig. 67 Fig. 74N, Fig. 79

Universal Pediatric Forceps

Root Pick Extracting Forceps



CSI-150 Fig.365/150SK Upper anteriors and Roots



CSI-151 Fig.360/17SK Lower Molars



CSI-152 Fig.370/151SK Lower anteriors and Roots



CSI-152 Forcep Fig.301



CSI-153 Forcep Fig.300

Extracting Forceps (American Pattern)

Oral Surgery

Perforated Light Weight Handles with Strong Serrated Gripping Action.




CSI-159
Fig. 1
Upper Anterior

Side View



CSI-163
Set of 10 Pieces
Fig. 1, Fig. 23, Fig. 74
Fig. 88L, Fig. 88R, Fig. 74
Fig. 106, Fig. 174, Fig. 859

CSI-154
Lower Molar
Fig. 23 Cow Horn




CSI-155
Fig. 35
Upper Premolar



CSI-156
Fig. 88L
Upper Molar Left




CSI-157
Fig. 88R
Upper Molar Right




CSI-158
Fig. 74
Lower Molars



CSI-160
Fig. 106
Lower Incisors




CSI-161
Fig. 174
Extra Narrow
Lower Anterior
Apical Retention
Forceps



CSI-162
Fig. 74
Upper Molars




Lower Roots




BABY EXTRACTING FORCEPS

Think Teeth

Keep up good dental habits at every age.



Pediatric Tooth Extracting Forceps Sets



CSI-164
English Pattern
Klein #137 Pedodontic,
Upper Incisors
Spring handle




CSI-165
English Pattern Klein #5
Pedodontic Lower Incisors.




CSI-166
English Pattern Klein #7
Pedodontic Lower Roots.




CSI-167
English Pattern Klein #6
Pedodontic Lower Molars,
Spring Handle



CSI-168
English Pattern Klein #139
Pedodontic Upper Premolars



CSI-169
English Pattern Klein #51S
Pedodontic Upper root




CSI-170
English Pattern Klein #3
Pedodontic Upper
Molars Universal



CSI-171
Set of 7 Pieces

Pediatric Tooth Extracting Forceps Sets



CSI-176
English Pattern
Klein #137 Pedodontic,
Upper Incisors
Spring handle



CSI-172
English Pattern Klein #5
Pedodontic Lower Incisors.




CSI-173
English Pattern Klein #7
Pedodontic Lower Roots.




CSI-174
English Pattern Klein #6
Pedodontic Lower Molars,
Spring Handle



CSI-175
English Pattern Klein #139
Pedodontic Upper Premolars




CSI-177
English Pattern Klein #51S
Pedodontic Upper root



CSI-178
English Pattern Klein #3
Pedodontic Upper
Molars Universal




CSI-179
Set of 7 Pieces



Extracting Forceps(American Pattern)

Oral Surgery



CSI-165
Fig. 62
incisors, premolars,
deciduous teeth



CSI-166
#10S
upper roots



CSI-167
Harris
Fig. 18R



CSI-168
Harris
Fig. 18L



CSI-169
Fig. 53R
upper molars, right




CSI-170
Fig. 53L
upper molars, left



CSI-171
NEVIUS
Fig. 88R
upper molars, right



08



CSI-172
Fig. 24
upper molars




CSI-163
210S
upper third molars




CSI-163
99C Kells
upper incisors,
canines, premolars




CSI-163
Nevius
Fig. 88L
upper molars, left



CSI-173
Tomes Fig. 69
upper and lower
roots




CSI-174
Hull Fig. 101
Deciduous teeth




CSI-175
Fig. 150 AS

Extracting Forceps(American Pattern)


Oral Surgery




09




CSI-180
Fig. 222
Lower wisdoms



CSI-176
Fig. 17
Lower molars




CSI-181
Fig. 286
Upper roots,
premolars and incisors



CSI-177
Fig. 16
Harris Lower molars,
COW HORN beak




CSI-182
Fig. 171
Wisdom and molar roots



CSI-178
Fig. 222
Lower wisdoms



CSI-183
Fig. 166
Universal Application




CSI-179
Fig. 69
Tones Fragments,
uppers and lowers




CSI-184
Fig. 169
Centrals laterals premolars
lower canines and all roots




CSI-186
Fig. 24
upper molars




CSI-185
Fig. 151
Cryer Lower premolars,
incisors, roots, universal



CSI-187
Fig. 76
Upper roots



CSI-188
Tweezer Albrecht



CSI-189
Scissor Handle
Lower Wisdoms


BONE RONGEUR

INSTRUMENTS FOR CUTTING BONES

- Most commonly used
- Sharp blades- squeezed together cutting/pinching through bone
- Leaf spring between the handle : instrument opens when hand pressure is released
- Repeated cuts without manually reopening




Oral Surgery



CSI-190
Blumenthal 90°
15cm



CSI-191
Blumenthal 45°
15cm



CSI-193
Blumenthal 30°
15cm



CSI-192
Friedman
14cm



CSI-194
Mini Friedman
14.5cm



CSI-195
Mini Friedman 90°
14.5cm

SOFT TISSUE NIPPERS



CSI-196 10cm
Compact, sharp nippers
with smooth action useful
for accentuating interproximal
gingival contour during
soft tissue surgery.

DOUBLE ACTION RONGEURS



CSI-197
Double Action
Side-Cutting
17 cm



CSI-198
Boehler
Double Action
STR

CSI-199 15 cm-2mm
CSI-200 15 cm-3mm
CVD

CSI-201 15 cm-2mm
CSI-202 15 cm-3mm



CSI-203
Beyer
Double Action
18 cm

KERISON RONGEUR




Quickly and Easily Enlarge
The Sinus Window For
Increased Access
& Visibility

CSI-204 -2mm
CSI-205 -3mm
Kerison Rongeur



HANDLES FOR ROOT ELEVATORS





Elevators Classification

- 1-According to use:
 - 1-Elevators designed to remove the entire tooth, straight elevators, hospital pattern and winter elevator
 - 2-Elevators designed to remove roots broken off at the gingival line e.g. Apexo elevator, Coupland and Ido lavien elevators
 - 3-Elevators designed to remove roots broken off half way to the apex e.g. curved elevator hospital pattern, winter elevator and Apexo elevator



HANDLES FOR ROOT ELEVATORS POINTS / A-B-C-D-E-F-G-H



11

Flexible Periotomes

Gold Titanium Coated Points With Micro Serration





CSI-206

Flexible Periotome 3mm STD








Root Screw Extractor




Flexitome (Set of 6 Pcs)

Super Flexibility, Micro Serrated, Special Alloy, Hard To Be Broken



Flexitome Cassette




Root Elevators

Oral Surgery


Apexo Root Elevators



Bein Root Elevators



Seldin Root Elevators




Apical Root Elevators



Cryer Root Elevators



Coupland Root Elevators



Crane Pick



Root Elevators

Oral Surgery

Lindo Root Elevators




Bernard Root Elevators



Claw-Lift Roots Elevator



Heidbrink Root Tip Picks



PDL Elevators

Oral Surgery

For those clinicians who are interested the entire luxating Elevators line.
has created sets with seven micro serrated instruments for atraumatic extraction.

PDL Elevators (Set of 7 Pcs)



PDL Elevators (Set of 7 Pcs)



PDL Elevators (Set of 7 Pcs)




Bernard Flexible (Set of 4 Pcs)



Root Elevators

Oral Surgery

Luxation Elevators




Luxation Elevators



Potts Luxation Elevators



Chompret Syndesmometer



Surgical Curettes



Lucas Surgical Curettes




Root Elevators

Oral Surgery

Laxation Set, Delicated Points



Ergonomic Handles Elevators Set




Dental Anesthesia Syringes

Oral Surgery

Self Aspirating Local Anesthesia Syringes



Suitable For Use with All 1.8ml
Dental Cartridges




CSI-340



CSI-341




CSI-342



CSI-343

Mechanical Aspirating Local Anesthesia Syringes


Suitable For Use with All 1.8ml
Dental Cartridges



CSI-344



CSI-345



CSI-346





CSI-347

Intraligamental Syringes



CSI-348



Made Of Stainless Steel Color Coated. Suitable For Use with All 1.8ml
Dental Cartridges. Inject 0.05ml with Every Click. Exact Dosage Every Time.
Supplied with Angled Tubes.

Gun Style Intraligamental Syringes

Use *Low's* Premium Quality Stainless Steel Syringe to give
the anesthesia to patients without any problem.



CSI-349

STRAIGHT
ANGLED


Made of Stainless Steel

Metric Thread
Inject 0.06ml with every click.
Exact dosage every time
Suitable for use with all 1.8ml dental Cartridges
Supplied with Straight and angled tubes

Handles for Mouth Mirrors

Oral Surgery

Handles for Mouth mirrors



18

Dental Mouth Mirrors



Double Sided Mouth Mirror

Front Surface



Matrix Holding Instruments




Dental Tweezers

- ▶ Tweezers are used particularly to take and place cotton pellets, wedges and articulation
- ▶ Tweezers are important Dental Lab Equipment, made of non-magnetic stainless steel.
- ▶ They are designed particularly to hold the tiniest article, provide better grip on the jaw, paper from the to the mouth and back.



Oral Surgery



19

Adson Tweezers (Tungsten Carbide)



Micro Tweezer



Lab Tweezers



DENTAL SCISSORS

Oral Surgery

Functions and Features:

- Scissors designed for suture removal.
- Blades – one has a curved edge and one has a notched edge.
- During suture removal, suture is scooped up by blade with the notched end before cutting to prevent tissue trauma.

General Scissors



Castroviejo Scissors


Extra Light Weight



CSI-416 15cm
CSI-417 17cm
Castroviejo Scissor, Straight

Made By Titanium

Lenox Flex-Spring Scissors



CSI-418 15cm
CSI-419 17cm
Castroviejo Scissor, Straight

CSI-420 Noyes 11cm Straight

CSI-421 Noyes 11cm Angular

Super Cut Scissors



CSI-422 Iris 11.5cm Straight Saw edge

CSI-424 Kelly 16cm Straight Saw edge

CSI-427 Beebee 11cm Saw edge, Straight

Scissors Tungsten Carbide

Tungsten Carbide scissors feature cutting edges of permanently welded Tungsten carbide inside. These Scissors will retain a fine sharp cutting edge for long time and will outlast any standard stainless steel dental scissors. Both handles are gold-plated for easy identification.



Needle Holders

Oral Surgery

Needle holders feature tungsten carbide inserts with a matched diamond-cut pyramidal profile for enhanced longevity, a secure, positive suture needle grasp and reduced risk of cutting the suture material. With a precision designed hinge, the smooth surface reduces the risk of catching and or cutting the suture material during the knot tying process, enhancing clinical efficiencies. For optimal grip and control, offers a variety of design options, including standard, finger and palm controlled closure.

Tungsten Carbide Needle Holders







Needle Holders



Titanium Made Instruments

Lightweight With A High Tensile Strength





A13-Single Ended

TITANIUM COATED
HANDLES (B)

HOLLOW HANDLES (A)



SOLID HANDLES (D)



SILICON HANDLES (E)



Explorers & Probes

Dental Diagnostic


Hollow Handls (A)



Titanium Coated Handls (B)




Solid Handls (D)



Silicon Handls (E)




Explorers Double Ended



Explorers Double Ended 2A


Probes Double Ended



Probe Double Ended 12


EXPROS

A double-ended instrument that combines two critical instruments required for every effective diagnostic set-up with an explorer on one end and a probe on the other, an Expros is excellent for dental examinations and charting.




Excavators

SPOONS: For Removal of Carious Dentin



Pluggers / Condensers

Used to Pack Filling Materials into Cavity Preparations



Hollow Handles (A)



Titanium Coated Handles (B)



Solid Handles (D)




Silicon Handles (E)



CARVERS

Used to Crave anatomical features and trim excess materials.




Burnishers / Hollenback / Carving Knives

Operative Instruments


Hollow Hands (A)



Titanium Coated Hands (B)



Solid Hands (D)




Silicon Hands (E)



Burnishers

Designed to condense, smooth, carve and polish amalgam.
Gold Titanium Coated Points.



PK-Thomas

Used for waxing procedures and techniques.



Carving Knives

Used for trimming excess filling material, flash and overhangs.



Restorative

Operative Instruments


Cord Packing Instruments

For atraumatic and accurate cord placement.




CSI-555 / BN1

Thin blade and rounded contour facilitates use in both thick and thin tissues without catching or dropping cord. Bilateral notch allows placement in limited access areas.




CSI-556 1 Yardley
Round non-serrated working end.



Serrated

CSI-557
CSI-1

Non-Serrated



CSI-558
7 Guyer




Serrated

CSI-559
113


Non-Serrated

Gingival Retractors


Protects tissue during cavity preparations and placement of composite restorations.




CSI-560 / 1 Meinershagen / For mandibular incisors




CSI-561 / 2 Meinershagen
For maxillary and mandibular premolars and canines.



CSI-562 / 3 Meinershagen
For maxillary central incisors and wide canines



CSI-563 / 4 Meinershagen
For all molars.




CSI-564
10 Goldman-Fox


Composite Filling Instruments

Also Available in Titanium Coated Colored Points.

Vertical and horizontal combination.



CSI-565



CSI-566




CSI-567




CSI-568



CSI-569



CSI-570



CSI-571



CSI-572




CSI-573




CSI-574
Woodson



CSI-575



CSI-576



CSI-577




CSI-578



CSI-579

Vertical in and offset angled.

Hollow Handles (A)



A1-Double Ended




A2-Double Ended



A3-Double Ended



A4-Double Ended



A5-Double Ended




A6-Double Ended



A7-Double Ended



A8-Single Ended



A9-Single Ended



A10-Double Ended



A11-Double Ended



A12-Double Ended



A13-Single Ended

Titanium Coated Handles (B)



B1-Double Ended



B2-Double Ended




B3-Single Ended




B4-Double Ended



B5-Double Ended



B6-Double Ended



B7-Single Ended




B8-Single Ended

Solid Handles (D)




D1-Double Ended




D2-Double Ended



D3-Double Ended



D4-Double Ended



D5-Single Ended



D6-Double Ended

Silicon Handles (E)



Restorative

Operative Instruments

Hollow Handles (A)



Titanium Coated Handles (B)




Solid Handles (D)




Silicon Handles (E)



Composite Filling Instruments



Composite Filling Instruments Set



Amalgam Carrier



Dental Scalers


Operative Instruments

Sickle Scalers

Designed for removal of supragingival calculus.


All Scalers & Cutters are Available in

J2-420 ECO	USA 440c Premium
Best Quality Medical Stainless Steel. High Chromium Content Good Strength Easy To Sharpe Economical Cost Available in Different Sizes.	440c High-end Stainless Steel High Chromium Content High Carbon Martensitic High Strength Better Edge Retention Titanium Coated Points




CSI-604

Contra angle design offers greater access to distal surfaces




CSI-605

Super-thin, contraangle design is for interproximal areas. Shank angulation provides accessibility to all anterior and premolar surfaces



CSI-607 / H5/33 - Hygienist-Jacquette


Super-thin sickle scaler combined with a straight sickle blade. For scaling interproximal and cervical areas of anterior teeth.



CSI-609


U15/33

Towner-Jacquette




CSI-611

For posterior supra-gingival interproximal scaling




CSI-613

Contra angle design offers greater access to distal surfaces




CSI-615

Nabraska 5/6




CSI-617

Nevi-1 / H5
Anterior



CSI-606


Hygienist-Towner



CSI-608


U15/30

Towner-Jacquette




CSI-610

204




CSI-612 Morse

Miniature blade for scaling crowded or overlapping anterior teeth excellent for stain removal




CSI-614

204S



CSI-616


NV1/H5



CSI-618

Nevi-2 / Posterior
A super thin posterior curved sickle scaler


Columbia Curettes



13/14 Columbia
(Regular)


Universal Application
CSI-619

13/14 Columbia
Rigid



CSI-620

2R/2L Columbia, Anterior application



CSI-621

4R/4L Columbia, Posterior application

Hollow Handles (A)



Titanium Coated Handles (B)



Solid Handles (D)



Silicon Handles (E)



Dental Currettes

Operative Instruments

Hollow Hands (A)



Titanium Coated Hands (B)




Solid Hands (D)




Silicon Hands (E)




Langer Currettes



Barnhart



Periodontal Files




Gracey Standard Curettes

Operative Instruments

- Area-specific to allow for deep scaling, root planing and periodontal debridement
- Offset blade provides a perfect working angulation for the tooth surface
- Offset blade angulation means only the lower cutting edge of each blade is used

Gracey Standard Curettes



Hollow Handles (A)



Titanium Coated Handles (B)



Solid Handles (D)



Silicon Handles (E)



Periodontal Knives

Operative Instruments


Hollow Handles (A)



Titanium Coated Handles (B)



Solid Handles (D)



Silicon Handles (E)




Periosteals

For reflecting and retracting the mucoperiosteum after incisions of the gingival tissue.




Goldman-Fox




CSI-659 Goldman-Fox

Spear Point – For cutting interdentally after initial incision. Severs buccolingually. Angulated blade for easier handling.




CSI-660 Goldman-Fox

Combination – Spear-pointed blade used in mandibular region from posterior direction. Modified sickle blade (elongated) used in interdental area for debridement.




CSI-661 Goldman-Fox

Spear Point – Similar to #6, except blades are offset for use in posterior region.



CSI-662

3/4 Buck




CSI-663

5/6 Buck

Dental General Instruments

Oral Surgery

Retractors & Mouth Gags



CSI-664
Infant, Small
Children, Medium
Adults, Large




Weitlaner
CSI-665 10.5cm
CSI-666 13.5cm




CSI-667
Stainless Steel
Double Ended Cheek
Retractors




CSI-668
Stainless Steel
Single Ended Cheek
Retractors




Cheek Retractors,
CSI-669 14cm
CSI-670 16 cm




Minnesota
CSI-671 Small
CSI-672 Large




CSI-674
Retractors
17 cm
Hooks 22 mm



CSI-675
Retractors, Cheek
retractors combination
21 cm, 2



LD17-114
140 mm




LD17-115
150 mm
CSI-676




SSC-664




CSI-677
Bishop



CSI-678
Austin Fig.1



Surgical Retractor
CSI-680 7x25mm
CSI-681 10x42mm
CSI-682 12x55mm




Implant
Retractor
23cm




Scalpel Handles




CSI-683
No. 5E




CSI-684
No. 5A




CSI-685
3K Mini



CSI-686
Fig.3



CSI-687
Fig.4




CSI-688
Fig.3G


Bone Files



CSI-689
Miller
Cross Cut




CSI-690
Miller 2X
Straight Cut




CSI-691
Wahl 1


Cement Spatulas & Wax Knives




CSI-692
Gritman




CSI-693
3



CSI-694
1 Medium
Cement



CSI-695
Large Wooden Handle
Wax Knives,




CSI-696
20 cm curved

Rubber Dam Instruments

Oral Surgery

Rubber Dam Forceps



















































01-Anisworth
02-Brewer
03-Clamps Stand
04-Rubber Dam Frame
05-Rubber Dam Clamp
Fig.212,211,8,8A
14,0,00,2,2A,7



01-Anisworth
02-Palmer
03-Clamps Stand
04-Rubber Dam Frame
05-Rubber Dam Guide
06-Rubber Dam Clamp
Fig.212,211,8,8A
14,0,00,2,2A,7



Rubber Dam Clamps

									
Fig. 15	Fig. 16	Fig. 24	Fig. 25	Fig. 26	Fig. 28	Fig. 30	Fig. 31	Fig. 51	Fig. 136
									
Fig. 139	Fig. 200	Fig. 201	Fig. 202	Fig. 203	Fig. 204	Fig. 205	Fig. 3	Fig. W3	Fig. 4
									
Fig. 56A	Fig. 7	Fig. 7A	Fig. 8	Fig. 8A	Fig. W8A	Fig. 10	Fig. 11	Fig. 12A	Fig. 13A
									
Fig. 14	Fig. 14A	Fig. 14AN	Fig. 22	Fig. 27	Fig. 207	Fig. 29	Fig. 206	Fig. 208	Fig. 209
									
Fig. 1	Fig. 1A	Fig. 2	Fig. 2A	Fig. W2A	Fig. 9	Fig. 00	Fig. 210	Fig. 211	Fig. 212
									
Fig. 6	Fig. 9	Fig. N06	Fig. U67	Fig. 214	Fig. P1	Fig. P2	Rubber Dam Clamp 9 pcs Kit Fig. 7,14A,2A,8A,4 13A,12A,9,2		

Dental Crown Instruments

Oral Surgery



CSI-702
Dental
crown remover
forceps




CSI-703
Crown Spreaders



CSI-704
WYNMAN
Crown Gripper, CARB-N-SERT




CSI-705
Key Crown Remover
Set of 3 pcs



CSI-706
Universal
Crown Remover
Hollow Handle



CSI-707
Automatic
Crown Remover
with 2 Tips




CSI-708
Crown Remover Set
Full Automatic



CSI-709
Boley Gauge




CSI-710
Napkin Holder
Silicon




CSI-711
Napkin Holder


Bone Gauges and Measuring Instruments



CSI-712



CSI-713




- It is The Most Frequently Used Bracket Positioning Aid.
- It is Made Of Stainless Steel.
- It Accurately Measures Height Of Bracket Placement From The Incisal Edge At 3.5, 4, 4.5 And 5mm.
- Flat Surface Rests On Incisal/occlusal Surface Of Tooth. Simplifies Seating The Bracket On The Teeth With Exactly The Right Distance Between Bracket Dot And Incisal Edge With Pencil Led Or Metal Tips.



CSI-714
Measuring Loop
Straight



CSI-715
Lwanson Crown
Gauge



35

Bone Pluggers



Sinus Lift Instruments




Bone Chisels



Tissue Punch Straight & Curved




Convexes Straight & Curved





CSI-732

Bone Ridge Mapping
Forceps



CSI-733

Measuring Caliper



CSI-734


Splinter Forceps
Straight



Splinter
Forceps
Curved

CSI-735


Bone Mill Plier



CSI-736


Molt 18.5cm
Auto Clavable
Plastic Tips.

Bone Scraper



CSI-737


Curved



CSI-738

Straight

Star Torque Wrench



CSI-739

Star Torque Wrench
6.35mm Universal



CSI-740

Bone Grafting Syringe
17 cm, Straight / Bend
I.D 2.5mm, O.D 3.5

Bone Grafting Syringe
17 cm, Straight / Bend
I.D 3.5mm, O.D 4.5

Bone Grafting Syringe
18 cm, Straight / Bend
I.D 4.5mm, O.D 5.5

Bone Grafting Syringe
18 cm, Straight / Bend
I.D 6.0mm, O.D 7.0



CSI-741

Suction
Tub



CSI-742

Fig. 6FR
Frazier
15 cm-Ø2mm



CSI-743

Debakery Atraumatic
Tissue Forceps Str Serrated
Teeth 16cm

SSC-713
Debakery
Atraumatic
Tissue Forceps
Cvd
Serrated Teeth
16cm



CSI-744


Bone Expander Kit

Bone Expander Kit
Expander: 2.6mm, 3.0mm
3.4mm, 3.8mm, 4.3mm
Ratchet Driver
Parallel Pin
Pilot Drill
Drill 1.4mm
Drill Driver 4mm
Hand Driver 4mm
Ratchet 4.0mm



CSI-745
Bone Crusher
Well

Periotomes




CSI-746

Micro Serrated

Dental Implants Instruments

Oral Surgery



CSI-747

Mini Drill Set of 6 Pcs.
in Stainless Steel Stand,




CSI-748

Standard Drill Set of 6 Pcs.
in Stainless Steel Stand,



CSI-749

Trephine Bur Set of 6 Pcs.
in Stainless Steel Stand,



CSI-750

Bone Compression
Prosthodontic Dental Kit
for Implant



CSI-751

Tissue Punch Kit
in Stainless Steel Stand



CSI-753

Trimmer Kit in
Stainless Steel Cassette,



CSI-754


Mini Expander Kit
In Stainless Steel Cassette



CSI-752

Bur Holder
for 12 Pcs

Bur Holder
for 06 Pcs




CSI-760

Conical Universal Kit
(Kit of 50 Pcs)
in Stainless Steel Cassette

Orthodontics Instruments


Oral Surgery

Distal End Cutters




CSI-762
Distal End Cutter
Universal cut
long Neck 12cm

Max: .021"(.54mm)
x .025"(.64mm)
Min: .012"(.31mm)




CSI-763
Distal End Cutter
Universal cut
long handles 14cm

Max: .021"(.54mm)
x .025"(.64mm)




CSI-764
Flush Cutter with
Hold 12cm

Max: .021"(.54mm)
x .028"(.71mm)
Min: .012"(.31mm)




CSI-765
Distal End Cutter
Universal cut
long handle 14.5cm

Max: .021"(.54mm)
x .025"(.64mm)
Min: .012"(.31mm)




CSI-766
Distal End Cutter
Universal cut
Box Joint 12cm

Max: .021"(.54mm)
x .025"(.64mm)
Min: .012"(.31mm)



CSI-767
Distal End Cutter
Universal cut
long neck 12cm


Max: .021"(.54mm)
x .025"(.64mm)
Min: .012"(.31mm)



CSI-768
Distal End Cutter
Universal cut heavy
pattern.Box Joint long
Handles 14cm


Max: .021"(.54mm)
x .025"(.64mm)
Min: .012"(.31mm)

Ligature & Pin Cutters




CSI-769
Mini Ligature & Pin
Cutters 15 degree
angled, Slim Beak
12.5cm

Max: .012"(.31mm)
Soft Wires




CSI-770
Ligature & Pin
Cutters wide beak
15 degree angled
12.5cm

Max: .012"(.31mm)
Soft Wires




CSI-771
Ligature & Pin
Cutters straight
12cm

Max: .012"(.31mm)
Soft Wires




CSI-772
Micro Ligature & Pin
Cutters 15 degree
angled 11.5cm

Max: .012"(.31mm)
Soft Wires




CSI-773
Ligature Cutters
45 degree angled
wide Handles
very Fine Tips.
11.5cm

Max: .012"(.31mm)
Min: .008"(.20mm)
Soft Wires



CSI-774
Ligature & Pin
Cutters, Box
Joint 1/3 Handles
Gold Platted 12.5cm


Max: .012"(.31mm)
Soft Wires




CSI-775
Micro Mini Ligature
Cutters 15 degree
angle 11cm

Max: .012"(.31mm)
Soft Wires


Pliers




CSI-777
Bracket Removing
Plier 12.5cm
straight




CSI-778
Band Removing
Pliers long neck
12.5cm



CSI-779
Slim Weingart
Plier 13cm




CSI-780
Weingart
Plier 12.5cm



CSI-781
Weingart
Plier 14cm





CSI-782
Nance Pliers
12.5cm




CSI-783
Four Step
Loop Forming
Pliers 13cm

Orthodontics Instruments


Oral Surgery



CSI-784	CSI-785	CSI-786	CSI-787	CSI-788	CSI-789	CSI-790	CSI-791	CSI-792
Three Jaw Pliers 12.5cm	Adrene Pliers 13cm	Bird Beak Pliers 12cm	Lab Bird Beak Pliers with serrated tips 14cm	Light Wire Bending Pliers long 12.5cm	Slim Bird Beak Pliers Tapered Beak 12.5cm	Rectangular Arch Forming pliers 12cm	Slim Arch Forming Pliers 12cm	Tweed Arch Forming Pliers 12cm
Max: .028" (.71mm) Wires	Max: .030" (.76mm) Wires	Max: .030" (.76mm) Wires	Max: .036" (.92mm) Hard Wires	Max: .020" (.51mm) Wires	Max: .022"x0.25" (.56mm x .64mm) Wires	Max: .022"x0.25" (.51mm x .64mm) Wires	Max: .022"x0.25" (.51mm x .64mm) Wires	Max: .022"x0.25" (.51mm x .64mm) Wires



CSI-793	CSI-794	CSI-795	CSI-796	CSI-797	CSI-798	CSI-799	CSI-800	CSI-801
Hollow chop Arch Forming Pliers 12cm	Hollow chop Arch Forming Pliers with 3 Grooves 12cm	Tweed Plier 12.5cm	Adhesive Removing Pliers	How Utility Pliers angled 12.5cm	How Utility Pliers straight 12.5cm	Tweed Pliers 12.5cm	Loop Forming Pliers with cutters. Tweed Type	Step Forming Pliers A-0.25 / B-0.50 C-0.75 / D-1mm
						Max: .022"x0.28" (.56mm x .71mm) Wires	Max: .025" Rectangular Hard wires	Max: .022" (.56mm) Wires



CSI-802	CSI-803	CSI-804	CSI-805	CSI-806	CSI-807	CSI-808	CSI-809	CSI-810
Ligature Forming Pliers	Niti Cinch Back Pliers, Slim Beak	Hook Cimping Pliers	Loop Forming Pliers, long Beak	Optical Pliers	Young Loop Bending Pliers	Crown Crimping Pliers	Johnson Contouring Pliers	Niti Cinch Back Pliers, 13.5cm Short Handle
	Max: .021"x0.25" (.54mm x .64mm) Wires		Max: .028" (.70mm) Hard Wires		Max: .028" (.70mm) Hard Wires			Max: .021"x0.25" (.54mm x .64mm) Wires

Clear Brace Aligner Pliers



CSI-811	CSI-812	CSI-813	CSI-814	CSI-815	CSI-816
Hole Punch Cutter	Tear Drop Cutter	Horizontal Dimple	Mini Dot Dimple	Vertical Dimple	Dimple Remover

Individual Torquing Pliers Set



CSI-817
Individual Torquing Pliers set
To make 6mm long torque 13cm
Can bend wire upto. .022" x .025"

Hard Wire Cutters



CSI-818
Hard wire Cutter



CSI-819
Hard wire Cutter




CSI-820
Hard wire Cutter



CSI-821
Hard wire Cutter


WIRE CUTTER TC GOLD

Tungsten Carbide (Tip Insert)




CSI-822 16 cm
CSI-823 18 cm
Wire Cutter
TC Gold
Double Action


Hand Instruments




CSI-824




CSI-825




CSI-826




CSI-827



CSI-828



CSI-829




CSI-830

Uses of Articulator

Primary purpose:

- 1) To hold opposing casts in a predetermined fixed relationship
- 2) To open and close.
- 3) To produce border and inter-border diagnostic sliding motions of the teeth similar to these in the mouth.



Technical Data:

Intercondylar distance 120 mm. / condylar inclination fixed at 30°
Bennet angle fixed at 7° / Incisal table inclination 10°
Distance between mounting plates 85 mm. / Weight 400 gr.
Opening 180°

PMC-887



Black

PMC-888



Blue

PMC-889



Yellow

Impression Trays

Oral Surgery

Impression Trays Dentulous Perforated



CSI-831

SIZE: S-M-L-XL

Impression Trays Edentulous Perforated



CSI-832

SIZE: S-M-L-XL

Impression Trays Dentulous Perforated



CSI-833

No Welding
Easy to Clean
Adjustable Front Wall

SIZE: S-M-L-XL

Impression Trays Edentulous Non-Perforated



CSI-834

SIZE: S-M-L-XL

Impression Trays Dentulous Perforated



CSI-835

No Welding
Easy to Clean
Adjustable Front Wall

SIZE: S-M-L

Impression Trays Edentulous Perforated



CSI-836

No Welding
Easy to Clean

SIZE: S-M-L

Impression Trays Dentulous Non-Perforated



CSI-837

SIZE: S-M-L

Lenox Impression Trays Sizes Chart

Uppers / Superior		Sizes / Measure (mm)				
		01 XS	02 XS	03 S	04 M	05 L
		56	60	66	71	73
		44	52	54	58	59
Lowers / Inferior		Sizes / Measure (mm)				
		L1 XS	L2 XS	L3 S	L4 M	L5 L
		56	66	70	73	76
		44	50	52	56	58

All Lenox Ins Impression Trays are Made by Surgical Grade Stainless Steel

Impression Tray Half Perforated



CSI-838

U1P / 12.5cm
For upper left jaw
For lower right jaw

U1P / 12.5cm
For upper right jaw
For lower left jaw



CSI-839

CSI-840

U1P / 12.5cm
Adjustable Perforated

U1P / 12.5cm
Adjustable Solid

Dental Hollow Wares

Oral Surgery




CSI-841
Cotton Holder



CSI-842
Spirit Lamp



CSI-843
Endodontic Box



Kidney Trays

With lid available	Without lid available
CSI-844 6"	CSI-848 6"
CSI-845 8"	CSI-849 8"
CSI-846 10"	CSI-850 10"
CSI-847 12"	CSI-851 12"




CSI-852 6" x 8"
CSI-853 8" x 10"
CSI-854 10" x 12"
With Lid



CSI-855 6" x 8"
CSI-856 8" x 10"
CSI-857 10" x 12"
Without Lid




CSI-858
Instruments Box
23x13x5cm




CSI-859
SCALER TRAY
10 x 20 mm




CSI-860
Dental PRF Box




CSI-861
Dental Scaler Clip Cassettes
for 10 Pcs




CSI-862 5 Pcs
CSI-863 B 7 Pcs
CSI-864 10 Pcs
CSI-865 20 Pcs



CSI-869
Dental Scaler Clip Cassettes
for 10 Pcs




CSI-870
Dental Scaler Cassettes
for 8 Pcs




CSI-866 5 Pcs
CSI-867 B 7 Pcs
CSI-868 C 10 Pcs



CSI-871
Dental Instruments Cassettes
for 5 Pcs



CSI-872
Dental Scaler Cassettes
for 10 Pcs



CSI-873 5 Pcs
CSI-874 7 Pcs
CSI-875 C 10 Pcs

Dental Scalers Cassettes, A, B, C, D

CSI-877
Dental Scalers Cassettes, A, B, C, D

Dental Instruments Cassettes

Oral Surgery



CSI-878

The Oral Surgery Kit 1001:

- 1: 3.0mm Aspiration Tip / 2: Serrated Lucas Surgical Curette, Small
- 3: Serrated Lucas Surgical Curette Large / 4: Miller #1 Bone File / 5: Cotton Pliers
- 6: Medium Straight Elevator / 7: #77R Back Action Elevator 8: Small Straight Elevator
- 9: Lower Universal Apical Forceps / 10: Lower Corshorn Forceps
- 11: Upper Universal Forceps / 12: Lower Universal Forceps / 13: Cone Socket Mirror Handle
- 14: Bard-Parker Scalpel Blade Handle / 15: Curved Kelly Haemostat / 16: #5 Cone Socket Mirror
- 17: Crile-Wood Needle Holder / 18: Seldin Flap Retractor / 19: Molt #9 Periosteal Elevator
- 20: UNC Perio-Probe / #23 Explorer / 21: Minnesota Retractor / 22: Blumenthal Rongeur
- 23: Super-Cut Kelly Scissors / 24: Adson/Brown Tissue Forceps, 2x8 Row Teeth
- 25: Extra Large Cassette



CSI-879

The Soft Tissue Kit # 2002

- Description
- 1: Perio Chisel #13 / 2: Large Cassette / 3: Locking Cotton Pliers
 - 4: Drop Blade Scalpel Handle / 5: Cone Socket Mirror Handle
 - 6: #5 Cone Socket Mirror / 7: Molt Curette 2x4 / 8: Crile-Wood Needle Holder / 9: Sugarman 3/4 Perio File / 10: Goldman-Fox #11 Perio Knife
 - 11: Perio-Probe / #23 Explorer
 - 12: Younger-Good 7/8 Curette
 - 13: Goldman Fox Scissors



CSI-880

Micro Surgery Kit# 3003

- 1: Perio Chisel #13L / 2: Perio Chisel / 3: Medium Cassette
- 4: Micro Cooley Forceps / 5: Twist Joint Castro-Viejo Needle Holder
- 6: Periosteal Elevator / 7: Castro-Viejo Scissors
- 8: Molt #9 Periosteal Elevator



CSI-881

The Implant Site Dev Kit # 4004

- 1: Small Distal Periosteal TT Coated / 2: Small Mesial Periosteal TT Coated
- 3: Serrated Lucas Surgical Curette, Small / 4: 3.5mm Graft Syringe
- 5: Large Cassette / 6: Upper Universal Deep Grip Forceps
- 7: Lower Universal Deep Grip Forceps / 8: Bone Condenser #12
- 9: Crile-Wood Needle Holder / 10: Mini Buser Periosteal Elevator
- 11: Super-Cut Iris Scissors / 12: Suture Forceps, 1.6mm



CSI-882

LD Hygiene Kit 5005:


- 1: Medium Cassette / 2: Cone Socket Mirror Handle
- 3: #5 Cone Socket Mirror / 4: UNC Perio-Probe, Color
- 5: #204S Scaler / 6: Gracey 11/12 Curette
- 7: Gracey 13/14 Curette / 8: Columbia 4R/4L Curette
- 9: Gracey 1/2 Curette / 10: H&S Scale



CSI-883

Basic Exam Kit 6006

- 1: Small Cassette / 2: Cotton Pliers
- 3: Cone Socket Mirror Handle
- 4: #5 Cone Socket Mirror
- 5: UNC Perio-Probe / #23 Explorer



CSI-884

Basic Exam Pack:7007

- 1: Regular #5 Mirror & Handle, Solid Handle
- 2: Explorer 17/23, Solid Handle
- 3: Mariani Cotton Forceps

INSTRUMENTS CARE

1STEP ONE:

Following clinical use, all instruments must be thoroughly rinsed with water to remove all residue prior to instrument processing.

STEP TWO:

Instruments should be placed into an ultrasonic bath for a period of 20 minutes to ensure complete instrument cleaning.

Note: Ultrasonic Solution should be changed daily or when the solution appears to be visibly soiled.

STEP THREE:

All instruments must be rinsed with water and completely dried before autoclaving. All stainless-steel instruments should be placed in separate sterilization pouches and not to be mixed with any carbide instruments during instrument processing.

NOTE: Only demineralized water should be used in autoclave equipment according to autoclave manufacturer instructions.

- ▶ All stainless-steel instruments must be cleaned and thoroughly dried prior to sterilization.
- ▶ We recommend using only a non-corrosive detergent.
- ▶ Abrasive brushes or materials should not be used to clean instruments.

2-STAINLESS STEEL INSTRUMENTS

We use many different types of stainless steels in the production of instruments dependent upon the design and application of the instrument. All stainless steels have a high nickel and chromium content to maximize corrosion resistance but will corrode and discolor when subjected to high concentrations of certain chemicals.

Do Not expose stainless steel dental instruments to the following chemicals. These chemicals will cause an adverse reaction and may destroy your instruments: Chlorine or Chlorinated products, Household Bleach, Tarter and Stain Remover, Aluminum Chloride, Aqua Regia, Barium Chloride, Bichloride of Mercury, Calcium Chloride, Carbolic Acid, Chlorinated Lime, Citric Acid, Dakin's Solution, Ferric Chloride, Ferrous Chloride, Hydrochloric Acid, Iodine, Lysol®, Mercury Chloride, Mercury Salts, Phenol, Potassium Permanganate, Potassium Thiocyanate, Sodium Hypochlorite (bleach), Stannous Chloride, Sulfuric Acid and Tartaric Acid (Tarter & Stain Remover)

Note: Be sure to only use solutions and chemicals that are compatible with stainless steel hand instruments.

Do Not use sodium chloride solution to clean the instruments. Never exceed temperatures 350o F / 177o C as this will have an adverse effect on the temper of the steel.

Recommended:

Steam sterilization in an autoclave is the benchmark process and is recommended. Sterilization requires a constant temperature of 134°C for 18 to 20 minutes.

3 TITANIUM NITRIDE COATED INSTRUMENTS

Titanium coatings are used to increase the surface hardness of instrument tips and reduce "pullback" when manipulating composite materials. Proper care should be taken to remove any residual composite materials from the blade within 5 minutes after use by wiping the blade with alcohol on a 2x2 gauze. Composite materials or residue can harden on the blade and affect the quality of future restorations. Composite residue left to harden on the blade can not be removed without damaging the surface finish of the composite placement instrument. Titanium Nitride Coated instruments can be cleaned and sterilized using the same recognized acceptable methods as stainless steel instruments

Note: Never use abrasives to clean the surface of any titanium coated composite placement instruments, elevators, periotomes or curettes.

HINGED INSTRUMENTS

All hinge instruments should be allowed to thoroughly dry in an open position after ultrasonic and manual cleaning procedures. We recommend that you only use lubricants such as surgical milk specifically formulated for dental and surgical instruments and follow manufacturer's instructions for applications.



COMFORT DENTAL INSTRUMENTS

Manufacturers & Exporters of Surgical & Dental instruments

☎ 646 688 5552

✉ info@cureventory.com

🌐 cureventory.com

📍 252 GlenHill Ave, Yonkes, NY 10701

